




APPLICATION NOTE

BRAINCUBATOR™


Extending the viability of acute brain slices

Yossi Buskila, Paul P. Breen, Jonathan Tapson, André van
Schaik, Matthew Barton & John W. Morley

Fresh brain slices have a limited lifespan, leading to experimental inconsistencies and unnecessary animal use. This study examines how a **controlled storage environment (BrainCubator) can extend slice viability for neuron isolation and electrophysiology recording**. By maintaining artificial cerebrospinal fluid (aCSF) at 16°C and minimizing bacterial growth through UVC sterilization, **stable recordings of pyramidal neuron resting membrane potential (RM) and action potentials (APs) were achieved for up to 36 hours**. Through the use of the Braincubator, slices retained physiological activity for longer, extending the reliable recording window and allowing researchers to maximize data collection from each preparation. By reducing extracellular solution temperature and inhibiting bacterial growth, the Braincubator helped preserve neuronal function, ensuring that layer V pyramidal neurons maintained stable resting membrane potentials and synaptic activity. This **extended viability reduced the need for additional animals, improving both research efficiency and ethical standards**.




SLICE VIABILITY – STANDARD CONTROL VS BRAINCUBATOR




Brain slice viability was assessed using propidium iodide (PI) and DAPI markers. Viability was measured as the ratio of live cells (PI negative) to total visible cells (DAPI positive + PI positive) in the rat cerebral cortex. The experiment involved imaging slices maintained either in the Braincubator or control conditions (n=54) at various time points (1- 36 hour intervals indicated above) after being sectioned at the same time. At 1 hour, no significant difference was observed in cell viability between conditions. However, **after 6 hours, slices maintained in the Braincubator showed significantly higher cell viability, with a maximal difference of 28% at 24 hours (45% vs. 17%, $p < 0.0001$).**

BACTERIAL CONTAMINATION

Bacterial numbers over time. Results display *Pseudomonas* and *Stenotrophomonas maltophilia* bacteria counts over time. These bacteria typically thrive in neutral pH and mesophilic conditions. In the Braincubator solution and slices (15–16°C, UVC filtration), bacterial growth was minimal compared to control conditions (22°C, no UVC). Bacterial growth was not observed in BrainCubator solution or slices until after 12 hours ($p < 0.0001$)





Electrophysiological properties of layer V pyramidal neurons were recorded in slices following incubation in the Braincubator. (a) I-V traces. Increasing step currents of 500 ms were injected into the soma through the recording electrode to reveal the input resistance and firing properties. (b) The resonance frequency was measured by injecting a chirp stimulation of 60 pA (peak to peak). (c) **Sample traces of spontaneous synaptic activity recorded from a neuron after 31 hrs in the braincubator, indicating extended viability of acute brain slices.**

	Control < 4 hrs (n = 10)	1–3 hrs (n = 12)	5–8 hrs (n = 7)	15–24 hrs (n = 9)	26–41 hrs (n = 9)
Resting membrane potential (mV)	-63 ± 2	-65 ± 2	-64 ± 1	-64 ± 1	-63 ± 1
Input resistance (MΩ)	176 ± 29	165 ± 18	163 ± 27	179 ± 21	145 ± 19
Membrane time constant (ms)	25 ± 4	28 ± 3	24 ± 7	25 ± 2	23 ± 3
Resonance frequency (Hz)	2 ± 0.3	2 ± 0.2	1.7 ± 0.3	1.5 ± 0.3	1.6 ± 0.3
First spike Amplitude (mV)	97 ± 5	97 ± 2	94 ± 4	93 ± 3	92 ± 3
Half-width spike amplitude (ms)	2.4 ± 0.2	2.5 ± 0.2	2.3 ± 0.1	2.6 ± 0.2	2.2 ± 0.1
mEPSC's frequency (Hz)	2.2 ± 0.3	1.75 ± 0.3	2.3 ± 0.4	1.5 ± 0.2	2.2 ± 0.3
mEPSC's amplitude (pA)	-7.3 ± 0.5	-6.4 ± 0.5	-7.8 ± 0.5	-5.7 ± 0.4	-8 ± 0.7

Electrophysiological properties of layer V pyramidal neurons were stable across incubation times in the Braincubator, with no significant differences between groups (one-way ANOVA). The control group (<3 hours) was included for comparison.

**Precisely regulates
pH, temperature,
and media**

**Tissue viability up to
36hrs**

**Halve the number of
animals sacrificed**

The Braincubator is an effective tool for significantly extending the viability of brain slices, **enabling neuron isolation and recording for up to 36 hours**. By combining reduced temperature (16°C) and UVC irradiation, it **minimizes bacterial growth and preserves slice integrity**, preventing neuronal damage typically caused by bacterial toxins. The use of the Braincubator enables **stable electrophysiological properties**, such as resting membrane potentials and functional synaptic activity, in neurons within the slice. This **controlled environment** not only **enhances the quality and duration of experimental data** but also **reduces the need for additional animal use**, making it a valuable tool in neuroscience research.

This study evaluated the effectiveness of the Braincubator™ system in maintaining acute brain slices over extended periods. Tissue viability, slice functionality, and bacterial contamination were assessed using fluorescence staining, electrophysiological recordings, and bacterial culturing. A standard interface chamber served as the control condition for comparison.

Animals

Male Wistar rats (2–5 weeks old) were housed in standard conditions with a 12:12 h light/dark cycle and ad libitum food and water. All procedures followed the University of Western Sydney Animal Research Committee guidelines (#A9452).

Slice Preparation and Incubation

Rats were anesthetized with 5% isoflurane, brains quickly removed and sectioned into 300 μm slices in ice-cold, carbogenated artificial cerebrospinal fluid (aCSF). Slices recovered for 30 minutes at 35°C in a submerged incubation chamber.

Experiment Setup: Slices were incubated in either the Braincubator™ (experimental) or a standard submerged chamber at room temperature (~22°C, control). with no active bacterial control or additional environmental regulation. Slices remained in their respective chambers for at least 30 minutes before further assessments.

Slice Viability and Electrophysiology

Slices were stained with propidium iodide (PI) and DAPI, and viability was assessed using confocal microscopy. Whole-cell patch-clamp recordings were conducted on layer 5 pyramidal neurons using a Multiclamp 700B amplifier. Resonance frequency was measured via impedance analysis using fast Fourier transform (FFT). Spontaneous miniature EPSCs were recorded in voltage-clamp mode at -70 mV in the presence of tetrodotoxin and picrotoxin.

Bacterial Detection

aCSF and brain slice samples were collected at multiple time points (1, 3, 6, 12, 24, and 36 h) from both incubation conditions and cultured on agar plates at 37°C. Colony-forming units (CFU) were counted, and bacterial identification was performed using MALDI-TOF mass spectrometry.




Braincubator™ System

The Braincubator™ – Recovery Incubation System

A custom incubation system was **designed to regulate pH, temperature, and bacterial contamination**. The system consisted of a main chamber housing brain slices and a separate UVC chamber for bacterial sterilization. The main chamber was insulated with polystyrene foam and monitored for pH and temperature. The UVC chamber contained a 5W Philips UVC lamp (254 nm), which operated intermittently (15–26 min every 15–30 min) to minimize heat buildup and prevent bacterial resistance. Solution was circulated between chambers via a peristaltic pump, and temperature was controlled using a Peltier thermoelectric cooler (TE Technology, MI).

Statistical Analysis

Data were analyzed using Student's t-tests or ANOVA, with results expressed as mean \pm S.E.M. Statistical significance was set at $p < 0.05$.

- 
- Buskila, Y. et al (2014): **Extending the viability of acute brain slices**. Scientific Reports 4, (2014).
- Cameron, M. et al (2016): **Calcium Imaging of AM Dyes Following Prolonged Incubation in Acute Neuronal Tissue**. PLoS One 11, e0155468 (2016).
- Jones, Roland S.G. et al (2016): **Human brain slices for epilepsy research: Pitfalls, solutions and future challenges**. Journal of Neuroscience Methods, Volume 260, 15 February 2016, Pages 221-232
- Cameron, M. A. et al (2017): **Prolonged Incubation of Acute Neuronal Tissue for Electrophysiology and Calcium-imaging**. Journal of Visualised Experiments, e55396 (2017).
- Cameron, M. A. et al (2017): **Differential effect of brief electrical stimulation on voltage-gated potassium channels**. Journal of Neurophysiology, doi: 10.1152/jn.00915.2016
- Sengupta, A. et al (2017): **Basolateral amygdala neurons maintain aversive emotional salience**. Journal of Neuroscience, 38:2460–17. doi: 10.1523/JNEUROSCI.2460-17.2017
- Bellot-Saez, A. et al (2018): **Astrocytic modulation of cortical oscillations**. Scientific Reports, 8:11565. doi: 10.1038/s41598-018-30003-w
- Gibson, Gabrielle D. et al (2018): **Distinct Accumbens Shell Output Pathways Promote versus Prevent Relapse to Alcohol Seeking**. Neuron, Volume 98, Issue 3, 2018, pp 512-520.e6
- Buskila, Y. et al (2019): **Extending the Life Span of Acute Neuronal Tissue for Imaging and Electrophysiological Studies**. Basic Neurobiology Techniques, pp 235-259 (Part of Neuromethods volume 152)
- Buskila, Y. et al (2019): **Dynamic interplay between H-current and M-current controls motoneuron hyperexcitability in amyotrophic lateral sclerosis**. Cell Death & Disease, 10, Article number: 310 (2019)
- Kekesi, O. et al (2019): **The differential impact of acute microglia activation on the excitability of cholinergic neurons in the mouse medial septum**. Brain Structure and Function, Volume 224, pp 2297–2309
- Kekesi, O., Buskila, Yossi (2020): **Method for prolonged incubation of brain slices**. Bio-Protocol, Volume 10, Issue 44 e3683

BRAINCUBATOR™

prolongs tissue slice lifespan, reduces animal usage, supports extended experimental time, and ensures more reliable experimental conditions

www.braincubator.com.au

