

made to measure

AUDIS-03/8M

Audio Monitor Module for the EPMS-07 System



Features:

- ⇒ Two operation modes
 - ⇒ PITCH mode: the input signal is converted into a tone
 - ⇒ NOISE mode: the input signal is highpass filtered to monitor the noise
- ⇒ THRESHOLD function sets the level for NOISE mode
- ⇒ THRESHOLD with 3 ranges (0.5V, 1V or 2V max.) and fine tuning in % of max. value
- ⇒ Signal input with BNC connector or via EPMS bus
- ⇒ Built-in speaker with volume control
- ⇒ HEAD PHONES 3.5 mm stereo jack audio connector



For more information contact:

General:

npi electronic GmbH
Phone: +49-7141-9730230
Fax: +49-7141-9730240
sales@npielectronic.com
www.npielectronic.com

North America:

ALA Scientific Instruments
Phone: +1-631-393-6401
Fax: +1-631-393-6407
sales@alascience.com
www.alascience.com

Great Britain:

Scientifica Limited
Phone: +44-1825-749933
Fax: +44-1825-749934
info@scientifica.uk.com
www.scientifica.uk.com

Switzerland:

Science Products Trading AG
Phone: +41-43-4880561
Fax: +41-43-4880562
info@science-products.com
www.science-products.ch

