

made to measure

ELP-01D

Electroporator for Tissue and Cultured Cells



Features:

- ⇒ Square wave pulse generation, isolated output
- ⇒ Single, multiple (up to 999) pulses, or continuous cycling operation
- ⇒ Digital cycle counter (000 to 999 cycles)
- ⇒ STOP and RESET function (switch operated)
- ⇒ DELAY and PULSE duration up to 9999.9 ms; resolution 0.1 ms
- ⇒ Manual, triggered (TRIGGER TTL) or gated (INPUT TTL) operation
- ⇒ AMPLITUDE setting with digital potentiometer (± 110 V, 1 V resolution)
- ⇒ Electrode resistance test (up to 9.99 k Ω , resolution 10 Ω)

For more information contact:

General:
npi electronic GmbH
Phone: +49-7141-9730230
Fax: +49-7141-9730240
sales@npielectronic.com
www.npielectronic.com

North America:
ALA Scientific Instruments
Phone: +1-631-393-6401
Fax: +1-631-393-6407
sales@alascience.com
www.alascience.com

Great Britain:
Scientifica Limited
Phone: +44-1825-749933
Fax: +44-1825-749934
info@scientifica.uk.com
www.scientifica.uk.com

Switzerland:
Science Products AG
Phone: +41-61-3033555
Fax: +41-61-3013570
info@science-products.com
www.science-products.com

