

made to measure

GND-BOX

Box for central setup grounding



Features:

- ⇒ 19 grounding points
- ⇒ 2 grounding points with cap nut
- ⇒ Internally connected by 1.5 mm, silver plated wires
- ⇒ Robust metal case (120 x 65 x 40 mm)

For more information contact:

General:

npi electronic GmbH

Phone: +49-7141-9730230

Fax: +49-7141-9730240

sales@npielectronic.com

www.npielectronic.com

North America:

ALA Scientific Instruments

Phone: +1-631-393-6401

Fax: +1-631-393-6407

sales@alascience.com

www.alascience.com

Great Britain:

Scientifica Limited

Phone: +44-1825-749933

Fax: +44-1825-749934

info@scientifica.uk.com

www.scientifica.uk.com

Switzerland:

Science Products Trading AG

Phone: +41-43-4880561

Fax: +41-43-4880562

info@science-products.com

www.science-products.ch