

*made to measure*


# ACI-01

for perfect chlorinating of silver wires



## Features:

- ⇒ Low-cost, fully automatic system for manufacturing Ag-AgCl electrodes
- ⇒ Makes the fabrication of chlorinated silver wires reliable and easy
- ⇒ Uses an external power supply
- ⇒ Works according to an optimized chlorinating procedure developed at the Max-Planck Institute in Göttingen, Germany
- ⇒ Automatic shut-off
- ⇒ Operation indicated by LED



ACI-01 working scheme

## For more information contact:

### General:

**npi electronic GmbH**

Phone: +49-7141-9730230

Fax: +49-7141-9730240

sales@npielectronic.com

www.npielectronic.com

### North America:

**ALA Scientific Instruments**

Phone: +1-631-393-6401

Fax: +1-631-393-6407

sales@alascience.com

www.alascience.com

### Great Britain:

**Scientifica Limited**

Phone: +44-1825-749933

Fax: +44-1825-749934

info@scientifica.uk.com

www.scientifica.uk.com

### Switzerland:

**Science Products Trading AG**

Phone: +41-43-4880561

Fax: +41-43-4880562

info@science-products.com

www.science-products.ch