

## Easy-to-use temperature controller



- Set temperature & GO
- Preset control parameters matched to heater
- Additional input for external sensor

### Choose one of these heaters:



Perfusion heaters



Heated microscope stage



Animal heating pads



Objective heater

### Why choose npi's temperature controllers:

- Designed for use in sensitive electrophysiological recordings
- **No switching artifacts** due to analog control
- **No electrical noise** due to linear power supply

#### Contact:

General:  
**npi electronic GmbH**  
Phone: +49-7141-9730230  
sales@npielectronic.com  
www.npielectronic.com

North America:  
**ALA Scientific Instruments**  
Phone: +1-631-393-6401  
sales@alascience.com  
www.alascience.com